

Wireless Device Alarm Management Solution IntelSTAT

Enhance Situational Awareness and Alarm Response Effectiveness

*One Platform... Multiple Integrated Applications
Customized, Scalable, Integrated Data Management in One Solution*

IntelNettReal-Time WirelessDevice Alarm ManagementCapability provides **-Flexible and Scalable Applications** - Customized to meet facility requirements for patient monitoring devices...

IntelSTAT - ImproveSituationalAwareness

Wireless Device Alarm Management

... for patient monitoring devices...



- Ventilator
- Pulse Oximeter
- IV Pump
- More...



Time	Room	Type	Resp Time	Unit	Status
08:00 PM	Room 101	Pulse/Ox	07:15	1	Control
08:15 PM	Room 102	Pulse/Ox	07:20	2	Control
08:30 AM	Room 103A	Pulse/Ox	08:45	3	Control
08:45 AM	Room 104	Pulse/Ox	08:50	4	Control
09:00 AM	Room 105A	Pulse/Ox	09:10	5	Control
09:15 AM	Room 106	Pulse/Ox	09:25	6	Control
09:30 AM	Room 107	Pulse/Ox	09:35	7	Control
09:45 AM	Room 108	Pulse/Ox	09:50	8	Control
10:00 AM	Room 109A	Pulse/Ox	10:10	9	Control
10:15 AM	Room 110	Pulse/Ox	10:20	10	Control
10:30 AM	Room 111	Pulse/Ox	10:35	11	Control
10:45 AM	Room 112	Pulse/Ox	10:40	12	Control
11:00 AM	Room 113A	Pulse/Ox	11:10	13	Control
11:15 AM	Room 114	Pulse/Ox	11:20	14	Control
11:30 AM	Room 115	Pulse/Ox	11:35	15	Control
11:45 AM	Room 116A	Pulse/Ox	11:50	16	Control
12:00 PM	Room 117	Pulse/Ox	12:10	17	Control
12:15 PM	Room 118	Pulse/Ox	12:20	18	Control
12:30 PM	Room 119A	Pulse/Ox	12:40	19	Control
12:45 PM	Room 120	Pulse/Ox	12:50	20	Control
01:00 PM	Room 121	Pulse/Ox	01:10	21	Control
01:15 PM	Room 122A	Pulse/Ox	01:20	22	Control
01:30 PM	Room 123	Pulse/Ox	01:35	23	Control
01:45 PM	Room 124	Pulse/Ox	01:50	24	Control
02:00 PM	Room 125A	Pulse/Ox	02:10	25	Control
02:15 PM	Room 126	Pulse/Ox	02:20	26	Control

• IntelSTAT can integrate to any device that has a nurse call port



IntelNettWireless Robust Capability

- Multiple audio-visual alarm awareness capability
- Scalable system (manages hundreds of devices)
- Wireless installation and infrastructure
- Alarms can be propagated through central station, wireless wallboards, pagers, radios, paging systems, Audio Visual Modules (AVM), Annunciation Panels and e-mail
- No signal interference with other systems

Manage and Document Events

- Detailed alarm status and history available upon demand
- Situational alarm response can be documented for staff oversight and as a risk management tool
- Data is stored in the system hard drive and in a cloud server
- Information can be accessed locally or remotely via email or other smart devices



IntelSTAT
Wireless Device Alarm Management Components

Alarm Communication Options

- Smartphone
- 2 Way Radio
- Digital Pager



Hallway Digital Marquee

- **Wireless Device Alarm Management for patient monitoring devices:**
 - Ventilator
 - Pulse Oximeter
 - IV Pump
 - and more...
- **We can integrate with any device that has a nurse call port**

Nursing Station

32" Dynamic Viewing Station



10" Remote Dynamic Viewing Console



IntelNett Data Management System (DMS)



Additional wireless system functionality:

- Fall Prevention sensor mats
- Wireless Incontinence Management System (WIMS)
- Resident Wandering

Hallway Digital Marquee



Nursing Office

10" Remote Dynamic Viewing Console



Seamless Email Reporting



Wireless Transmitter & Device



Room 101

Wireless Transmitter & Device



Room 103

Wireless Transmitter & Device



Room 105

Wireless Transmitter & Device



Room 107



PLUS! Upgrade Your IntelNett System to include:

- **Wireless Fall Protection Solutions**
- **Wireless Incontinence Management System (WIMS)**
- **Wireless Patient Wandering and Asset Tracking Solutions**
- **Wireless Nurse Call and Mobile Nurse Call Solutions**

